
Harmonic™



Ethicon Generator G11

For HARMONIC™. For ENSEAL™. For the future.

Ethicon Generator G11:

For HARMONIC™. For ENSEAL™. For the future.

The Ethicon Generator G11 combines advanced technology, multifunctionality and touch screen simplicity in one sleek and easy-to-operate unit. No matter which Ethicon advanced Energy device you plan to use, the Ethicon Generator G11 will be the only power source you'll need.



- Intuitive Touch Screen**
- for fast and easy setup

Example of ENSEAL™ screen



Future Proof

- Upgradeable software via USB memory stick



Small and compact

- to take up less space in today's operating rooms

Universal Connector for Ethicon Energy Ultrasonic and Advanced Bipolar Platforms

- with automatic instrument recognition

Simple Design. Simple Setup*.

The shape of the hand-piece only permits instruments to be inserted correctly. The magnetic assist ensures that instruments are fully seeded into the receptacle. Whether you plug in a HARMONIC™ or ENSEAL™ device the generator will automatically go to the correct settings.



Specialties:

HARMONIC™ and ENSEAL™ technology is used in a variety of surgical specialties including: General, Gynecology, Colon/Rectal, Bariatric, Thoracic, Head & Neck, Urology, Orthopedics and Plastics. Ethicon is a proven leader in advanced energy with greater than 33 million procedures worldwide.^{1**}



Ordering:

Code	Description	Quantity/Sales Unit	Compatible with
GEN11	Ethicon Generator G11	1	All HARMONIC™ and ENSEAL™ devices
FSW11	Generator Foot Switch	1	Ethicon Generator G11
CRT11	Generator Cart	1	Ethicon Generator G11
HGA11	HARMONIC™ Connector	1	Ethicon Generator G11 and some HARMONIC™ hand pieces (HP054 & HPBLUE)
EGA11	ENSEAL™ Connector	1	Ethicon Generator G11 and current ENSEAL™ disposable instruments*

* except NSLG shears that do not require the ENSEAL™ Connector

** Internal global sales data as of 30 Sept 2020

1. Ethicon, 30092020, Energy Historical Units, Sept 2020, Data on File (123485-211124)

Please always refer to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

Ethicon Endo-Surgery (Europe) GmbH
Hummelsbütteler Steindamm 71
22851 Norderstedt, Germany

www.jnjmedtech.com

ETHICON

Johnson & Johnson SURGICAL TECHNOLOGIES

© Ethicon Endo-Surgery (Europe) GmbH 2022, 055117-220923 EMEA