

Innovative Usability for Enhanced Safety



Drive Gas Free

With electronically driven by high-performance turbine



8-Inch Capacitive Touch Screen

For intuitive operation and efficient interaction



Automatic Leak Test

Graphical guidance and one-click start



Inspiration Hold Function

Prevents respiratory motion interference with CT imaging and surgery

Maintain the plateau pressure

Prevent respiratory motion interference

Configurable start delay & hold duration



Veterinary Anesthesia Super Hub



* Configuration details subject to actual order specifications



LinkedIn Facebook Instagram

Q | Search Mindray Animal Medical

www.mindrayanimal.com
service@mindrayanimal.com

P/N:ENG-Veta 5Plus-210285X6P-20240802
 ©2024 Shenzhen Mindray Animal Medical Technology Co., Ltd. All rights reserved.

Veta 5 Plus Veterinary Anesthesia Machine

Redefining Veterinary Anesthesia Machines





Empowering Safety Elevating Confidence

Since its inception, the Veta series of veterinary anesthesia machines has consistently provided safer and more user-friendly anesthesia methods for clinical use.

The Veta 5Plus integrates the concept of active safety into veterinary anesthesia, innovatively developing the high-level ventilation mode—VS+*, which combines the experience of clinical experts with Veta's well-tested ventilation strategies to offer clinicians an intelligent anesthesia experience comparable to L3 assisted driving.

Innovative VS+ A New Era of Active Safety

Active Ventilation Adjustment

Corrects EtCO₂ by adjusting mechanical ventilation parameters



Instant Response

Synchronized adjustments guided by real-time data



Target EtCO₂ Control

Continuously tracks patient EtCO₂



VS+



Classic VS Superior Safety in Ventilation



Breathing Synchronization

Intelligent trigger for reduced respiratory resistance



Backup Ventilation

Effective prevention of hypoventilation with automatic backup ventilation

VS



Smart Pre-set

Automatic adjustment of mechanical ventilation parameters based on patient weight



Data Interconnection



Boundary-Breaking Connectivity Anesthesia Super Hub

Interconnecting data between monitors and anesthesia machines, The Veta 5Plus directly reads the patient's EtCO₂ from monitor and applies the VS+ ventilation strategy

Ensuring Active Safety

Based on VS Mode

Developed from well-tested VS patterns

Multiple Safeguards

Adjusts within safety ranges for volume/pressure/frequency

Gradual Adjustment

Rapidly adjusts when EtCO₂ exceeds the set range

Highest-Level Alarms

Activates if EtCO₂ is not corrected, alerting for immediate intervention



*Before using VS+, please thoroughly read the user manual and continuously monitor the depth of anesthesia and the patient's condition during use.