



OBSERVA SERIES

PrimaLED

The first and only examination lamp with the possibility to choose an ambient light in addition to traditional concentrated light. The Ergo-Spring balancing system makes PRIMALED very easy to handle and stable. PRIMALED is ideal for any type of installation, from the outpatient department to the intensive care unit.

Fixing Systems





PRIMAFLEX

Gooseneck Arm

PRIMA

Joint Arm

Performances

Primaled

Light intensity at 0,50 m distance (Ec)	105.000 lux
Light head diameter	19,5 cm
Color temperature	2 selections: 4.000 - 4.500K
Color rendering Index (CRI)	95 Ra
d10 light field diameter where illuminance reached 10% of Ec	150 mm
Average LED life	> 60.000
Control of the illuminance	25 - 100 %
Electrical absorption	10 W - 20 VA

Directive 93/42/EEC (main directive and further amendments) - Norm IEC 60601-2-41